



Product Service

CERTIFICATE

No. Z1 042097 0071 Rev. 00

Holder of Certificate: **Invacare Rea AB**
Org. No. 556098-9211
 Växjövägen 303
 343 71 Diö
 SWEDEN

Certification Mark:



Product: **Medical Beds**
Medical Bed for children

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition, the certification holder must not transfer the certificate to third parties. This certificate is valid until the listed date, unless it is cancelled earlier. All applicable requirements of the testing and certification regulations of TÜV SÜD Group have to be complied. For details see: www.tuvsud.com/ps-cert

Test report no.: 713184357

Valid until: 2026-01-31

Date, 2021-02-25

(Torsten Zimmer)



Product Service

CERTIFICATE

No. Z1 042097 0071 Rev. 00

Model(s): NordBed Kid Small
NordBed Kid Medium
NordBed Kid Large

Brand Name: INVACARE

Parameters:

Supply voltage	100-240 VAC
Frequency	50/60 Hz
Maximum current	2.5 A
Mode of operation	10% (2 min. on / 18 min. off)
Protection class	II
Applied parts	Type B
Degree of protection	IPX6
Application environment	3 & 4
Max. patient weight	70 kg
Safety working load	small: 100 kg, medium: 115 kg, large: 115 kg
Width	90 cm
Length	small: 146 cm, medium: 166 cm, large: 185 cm
Range of the height of the mattress support platform	30 – 80 cm

Tested according to: IEC 60601-1:2005
IEC 60601-1:2005/AMD1:2012
IEC 60601-1-6:2010
IEC 60601-1-6:2010/AMD1:2013
IEC 62366-1:2015
IEC 60601-1-11:2015
IEC 60601-1-2:2014
EN 50637:2017

ZERTIFIKAT ◆ CERTIFICATE ◆ CERTIFICADO ◆ CERTIFICAT ◆ 認證證書 ◆